



A Clinician's Guide to Restoring Screw Retained Full Zirconia Implant Crowns



www.AbsoluteDentalLab.com

NobelProcera® Full Contour Zirconia (FCZ)

What is NobelProcera® FCZ?

This innovative technology allows fabrication of a monolithic zirconia screw retained implant crown. The process eliminates the need for cement and therefore eradicates possible cement residue. This has been shown to be a common cause of inflammatory soft tissue complications* and in most cases the number one reason for implant failure. The system features a metal adapter into the implant creating a very strong and lasting connection.



FCZ System Features

- Extremely efficient and stress-free chairside protocol
- No cementation equates to no implant failures due to cement residue
- Available on all Nobel Conical Connection implants
- Design file archiving allows for impression-less remake if ever required
- Metal Adapter creates strong connection between implant and abutment-crown
- Monolithic design minimizes chipping under function
- Compatible with all Vita shades


Request Nobel Conical Connection implants to allow for this technology!

NobelProcera® Angulated Screw Channel Correction (ASC)

What is NobelProcera® ASC?

This technology enables the restorative team to re-angulate the screw access channel without the use of a correction abutment. Combining the ASC & FCZ technologies allows us to design Hybrid Anterior as well as Monolithic Posterior zirconia screw retained crowns. ASC allows the surgeon to place the implant in a more ideal position without having to adjust for access. This allows for a screw retained option without compromise.

ASC technology allows us to
re-angulate screw access
without the use of a
correctional abutment



FCZ System Features

- Allows for a more ideal implant placement without the need to consider screw access
- Enables screw retained options in both anterior and posterior applications
- Adjust screw access up to 25 degrees
- Rotate access up to 360 degrees to manage accessibility
- Traditional Nobel Biocare torque settings
- Omnigrip™ screwdriver fits standard Nobel Biocare torque wrench

Full Contour Zirconia Case Flow


Case courtesy of Dr. Thomas Leech

Posterior Screw Retained
Monolithic Crown


IMPRESSION AND PROCESSING



Implant level impression



Absolute digital design



NobelProcera® CAD milling & sintering



Absolute Lab final contouring & staining

Available on all Nobel Biocare Conical Connections

DELIVERY



Attach metal adapter



Position the crown



Adjust interproximal & occlusal contacts



Torque to Nobel Biocare protocols



Protect screw head with cotton plug



Plug access using composite material



Final delivery

"I delivered my first "FCZ" implant retained crown in less than 10 minutes. Prescribing this procedure with its retrievable option and no cement residue will change how we restore implants in the future. I truly believe this technology to be a more predictable, economical and functional way to restore dental implants."

– Thomas A Leech DDS, Dentistry at the Park, NC

Combination Full Contour Zirconia with ASC

Case courtesy of Dr. Brandon Kofford

Hybrid Monolithic/Layered
option for Anterior and
Posterior application


IMPRESSION AND LAB PROCESS



Removal of Immediate
Temporary




Implant level impression



Diagnostic wax-up



Digitizing the wax-up on the
NobelProcera® 2G Scanner



Digital cutback and access
correction using ASC



NobelProcera® Milling
and Sintering



Cut back and preparation for
layering Ceramics



Absolute Lab Processing -
Porcelain application to facial



Lab delivery of
processed crown

DELIVERY APPOINTMENT



Position and finger
tighten crown



Radiographical verification of
seating



Test and adjust occlusal and
interproximal contacts



Torque to traditional Nobel
Biocare torque settings



If using ASC a special screw will
be provided with case. Note
the blue screw head for ASC.



A special Omnigrip™ screw
driver is needed to engage an
ASC screw. This driver fits into a
standard Nobel Biocare wrench.



Protect screw head



Close access with composite



Monolithic Lingual /Layered Facial



Final Delivery

"After restoring my first hybrid FCZ screw-retained crown with Angulated Screw Channel Correction, I found the process very predictable and the results achieved to be both functional and esthetic. By prescribing these technologies I was able to solve the cement residue and screw access issues in one application."

- Brandon Kofford DMD, MS, FACP, Royal Oak Dental Group, NC

Clinical & Laboratory Components needed to restore FCZ & ASC

Noble Biocare order numbers shown, contact your rep for more information.

Impression Coping*

NP: 36258 narrow x short
36260 narrow x long
36259 intermediate x short
36261 intermediate x long

RP: 36263 narrow x short
36262 narrow x long
36265 intermediate x short
36264 intermediate x long
36267 wide x short
36266 wide x long

WP: 37855 intermediate x short
37854 intermediate x long
37857 wide x short
37856 wide x long

CLINICAL NEEDS

Unigrip – manual driver

29148 20mm // 29149 28mm // 29150 36mm

Omnigrip – manual driver – ASC application

37367 20mm // 37377 28mm // 37378 36mm

ASC Driver for torque wrench

37379 20mm // 37380 25mm // 37381 30mm // 37382 35mm

Nobel Biocare torque wrench

34584

Prosthetic Kit

3744

Lab analogs

NP: 36697 // RP: 36698 // WP: 37879

FCZ wax-up sleeve

NP: 37449 // RP: 37450 // WP: 37608

Abutment wax up holders

NP: 36746
RP: 36747
WP: 37568

NobelProcera® 2G scanner

37260

LAB NEEDS



**Narrow, intermediate or wide refers to emergence profile.
Long (14mm) or short (10mm) refers to length of post.*

Absolute Dental Lab

Established in 1994, Absolute Dental started as a fixed prosthetics lab serving clinicians in the Triangle area of North Carolina. Two decades later, Absolute's restorative focus is much broader but their attention to product detail and exceptional customer service has not changed.

Today, the Absolute team is renowned for their expertise in creating world-class dental esthetics. Their use of cutting edge technology in CAD and milling departments, as well as their extraordinary dental implant and high end removables sections, enables them to deliver lifelike and functional dental prosthetics.

Staying abreast of new technologies, yet only implementing relevant protocols and procedures, has earned Absolute a reputation for being a trusted partner to discerning clinicians throughout the United States.

Serving their customers with Absolute Excellence has always been the primary focus of the owners, branch partners and team members...this remains true today.



Conrad J. Rensburg

ND & NHD in Tech.

Conrad J. Rensburg graduated under full scholarship from Pretoria Tech in 1992, with a four-year Baccalaureate Degree in Technology. Since that time, he has specialized in fixed prosthetics with a heavy emphasis on dental implant restorations.



While employed by ADW Dental as Senior Ceramist, Mr. Rensburg received managerial certification from the SADTC in 1995. After, he owned and operated RensTech Fixed Prosthetics in Pretoria, South Africa. In 1999 he relocated to the U.S. and was hired as the General Manager of Absolute Dental Services. During that time, he created the nationally renowned Absolute Implant Department.

Mr. Rensburg and his business partner, Drew Van Aarde, purchased Absolute Dental in 2004. With Drew as Senior Ceramist and Conrad as head of Absolute Dental Implant and Pre-Ceramics sections, they created a world class dental lab in the Triangle area of North Carolina.

Mr. Rensburg is a board-certified technician in good standing with the SADTC and NCDLA. He has been a CE accredited speaker for several dental implant manufacturers since 2002, and has been the keynote speaker at special events across the U.S. His seminars are focused on educating restorative dentists on the latest techniques and materials, as well as keeping them up-to-date with the ever-evolving dental implant market.



www.AbsoluteDentalLab.com