



ABSOLUTE DIGITAL DENTURES



IN-HOUSE TURNAROUND

STANDARD DELIVERY

- 2 Days – Custom Tray
- 2 Days – Bite Rim
- 3 Days – Printed Tooth Try-In
- 3 Days – Final Delivery

EXPEDITED DELIVERY

- 24-hour turnaround from tooth try-in to final denture delivery

Must be requested through email. No phone call requests will be honored
Request expedited case through: expedite@absolutedentalservices.com
Only available if no changes to printed tooth try-in are required
Delivered within 24 hours from receipt of email requesting expedited processing
Not available in all areas – call for info before scheduling



CLINICAL PROCEDURES

STANDARD CLINICAL PROCEDURE

- 1st appointment: Initial impression. 2nd: Custom tray impression. 3rd: Bite registration.
- 4th: Printed tooth try in. 5th: Final delivery

DIGITAL SCAN OR DICOM

- 1st appointment: Hard reline denture. Scan with CBCT or digital scanner. Send STL, opposing and bite.
- 2nd: Tooth try-in. 3rd: Final delivery

EXISTING DENTURE

- 1st appointment: Hard reline existing denture. Send denture, opposing & bite - lab will return existing denture within 24 hours. 2nd: Tooth try-in or final delivery

IMMEDIATE DENTURE

- 1st appointment: PVS upper and lower impression with bite. 2nd: Final delivery
- 1st appointment: Digital impression. 2nd: Final delivery

PRODUCT INFORMATION

All digital dentures carry a one time, 12-month replacement warranty (see warranty details). All categories of Absolute digital dentures are processed using only premium multi-layer teeth. Denture bases are processed on state of the art, 3D Carbon printers and utilize Dentsply high impact, high strength Lucitone 3D printable base material. Pink tissue is available in five shades. Premium teeth are fused into the denture base using a proprietary FDA cleared cohesion process.